

DIAGNOSTICO DE PLEURITIS TUBERCULOSA CON ADA

| | Diagnostico TBC (+) | Diagnostico TBC (-) | Total |
|--------------|---------------------|---------------------|-------|
| Test ADA (+) | 17 | 2 | 19 |
| Test ADA (-) | 28 | 53 | 81 |
| Total | 45 | 55 | 100 |

-Sensibilidad: $a / (a+c) = 17 / (17+28) = 0.377 = 37\%$

-Especificidad: $d / (b+d) = 53 / (2+53) = 0.963 = 96\%$

-VPP: $a / (a+b) = 17 / (17+2) = 0.894 = 89\%$

-VPN: $d / (c+d) = 53 / (28+53) = 0.654 = 65\%$

-Sensibilidad/ 1- especificidad= $0.37/1-0.96 = 9.25$.